

Model: CMF86V80

-86°C Ultra Low Temperature Freezer

Overview

- -40°C to -86°C temperature range
- External digital controller and temperature display
- Annual service & Calibration available
- Fitted lock

2 YEARS PARTS & LABOUR WARRANTY

ANNUAL SERVICE & CALIBRATION AVAILABLE



Key Features



Controller

Controller system. Read, record, reset with the touch of one button. Adjustable temperature setting.



Flexible storage

Multiple shelf locations. Back stop on every shelf provides a flexible storage solution.



Spark-free

No internal light or fan, suitable for a spark free laboratory.

Details

Model: CMF86V80

Gross Capacity (litres)	86
Net Capacity (litres)	80
Unit External Width (mm)	665
Unit External Depth (mm)	720
Unit External Height (mm)	1265
Product Weight (kg)	85
Packed Width (mm)	685
Packed Depth (mm)	730
Packed Height (mm)	1335
Packed Weight (Kg)	95
Temperature Range	-40°C to -86°C
External Finish, Internal Finish	Painted Steel, Stainless Steel
Battery Back Up	No
Door Finish, Door Type, Door Style	Painted Steel, Solid, Hinged, Right Hand
Open Door Sensor	No
Lock	Yes
Adjustable Feet (number)	Yes, 2
Porthole	No
Porthole Location	N/A
Temperature Display on Panel, Temp Alarm, Remote Alarm Contacts, Audible Alarms, Visual Alarms	Yes

Number of Shelves	1
Number of Shelf Positions	6+
Shelf Material	Stainless Steel
Interior Light	No
Light Type	N/A
Light Activation	N/A
Fan	No
Spark Free	Yes
Number of Probes	1
Probe Type	Air
Ambient Temperature Range	+10°C to +30°C
Refrigerment Type	Mixed
CFC Free, HCFC Free, Ammonia Free	Yes
Voltage Supply	220-240V/ 50Hz
Certification	CE/UK
Noise Level (dbA)	56
Input Power (W)	450
Cable Length (Cm)	20
Energy Consumption (KWh/annum)	7.7
Energy Rating	F

Related Products



CMLFR145
145L Spark Free
Lab Fridge



CMLFZ115
115L Spark Free
Lab Freezer



CMLFF260
260L Combi
Fridge/Freezer