

Model: CMLFF260

## 260L Combi Fridge/Freezer

### Overview

- +2°C to +10°C & -16°C to -24°C temperature range
- Core temperature display, represents your vaccine vial through **iCOOL** temperature probe (fridge only)
- Freezer has an air probe
- External digital controller and temperature display
- Battery backup on both Fridge and Freezer sections
- Annual service & Calibration available
- Fitted lock with 2 keys



**5 YEARS PARTS & LABOUR WARRANTY**


**ANNUAL SERVICE & CALIBRATION AVAILABLE**

### Key Features



#### Controller

Controller system. Read, record, reset with the touch of one button. Adjustable temperature setting.



#### Porthole

Factory fitted porthole enables quick and easy access for any secondary temperature monitoring probes.



#### Battery

Backup should there ever be a power outage the battery backup will continue to monitor the internal temperature of the freezer.



#### Spark-free

No internal light or fan, suitable for a spark free laboratory.



#### Flexible storage

Multiple shelf locations. Back stop on every shelf provides a flexible storage solution.

## Details

Model: CMLFF260

<b>Gross Capacity (litres)</b>	150 Fridge & 130 Freezer
<b>Net Capacity (litres)</b>	145 Fridge & 115 Freezer
<b>Unit External Width (mm)</b>	595
<b>Unit External Depth (mm)</b>	595
<b>Unit External Height (mm)</b>	1645
<b>Product Weight (kg)</b>	79.5
<b>Packed Width (mm)</b>	650
<b>Packed Depth (mm)</b>	650
<b>Packed Height (mm)</b>	1740
<b>Packed Weight (Kg)</b>	85.5
<b>Temperature Range</b>	+2°C to +10°C -16°C to -24°C
<b>External Finish, Internal Finish</b>	Painted Steel, Painted Steel
<b>Battery Back Up</b>	Yes on Both
<b>Door Finish, Door Type, Door Style</b>	Painted Steel, Solid, Hinged, Right Hand
<b>Open Door Sensor</b>	Yes on Both
<b>Lock</b>	Yes on Both
<b>Adjustable Feet (number)</b>	Yes, 2
<b>Porthole</b>	Yes on Both
<b>Porthole Location</b>	Rear, Middle
<b>Temperature Display on Panel, Temp Alarm, Remote Alarm Contacts, Audible Alarms, Visual Alarms</b>	Yes

<b>Number of Shelves</b>	4 Fridge, 3 Freezer Baskets
<b>Number of Shelf Positions</b>	13
<b>Shelf Material</b>	Plastic Coated Wire
<b>Interior Light</b>	No
<b>Light Type</b>	N/A
<b>Light Activation</b>	N/A
<b>Fan</b>	No
<b>Spark Free</b>	Yes
<b>Controller Type</b>	Digital, External
<b>Number of Probes</b>	3 (2 in fridge, 1 in freezer)
<b>Probe Type</b>	Air & Load (Fridge), Air (Freezer)
<b>Ambient Temperature Range</b>	+10°C to +30°C
<b>Refrigerant Type</b>	R600a
<b>CFC Free, HCFC Free, Ammonia Free</b>	Yes
<b>Voltage Supply</b>	220-240V/ 50Hz
<b>Certification</b>	CE/UK
<b>Noise Level (dba)</b>	65
<b>Input Power (W)</b>	120
<b>Cable Length (Cm)</b>	16
<b>Energy Consumption (KWh/annum)</b>	1.84

## Related Products



**CMLFZ115**  
115L Spark Free  
Lab Freezer



**CMLFZ47**  
47L Spark Free  
Lab Freezer



**CMF86V80**  
-86°C Ultra Low  
Temperature Freezer

[www.coolmed.co.uk](http://www.coolmed.co.uk)

0161 772 5666 | [info@coolmed.co.uk](mailto:info@coolmed.co.uk)

Unit 8 Fourways Trading Estate, Trafford Park, Manchester, M17 1SW

**COOLMED**  
coolmed.co.uk